

Hill-Rom Overbed Table



At a glance

- Four different models provide a perfect solution for every care environment
- Optional food tray makes meal time convenient
- Three base options provide compatibility with every Hill-Rom bed



PatientMate® Overbed Table with storage compartment
Available with Standard or Low Base (OBT230)



PatientMate® Jr. Overbed Table
Available with Standard or Low Base (OBT220)




Standard Overbed Table
Available with Standard or Low Base (OBT210)



C-base Overbed Table (OBT100)


Hill-Rom Overbed Table



PatientMate® Overbed Table



PatientMate® Jr. Overbed Table



Standard Overbed Table



Standard Base

Low Base

Choose from two base heights for PatientMate®, PatientMate® Jr. and Standard Overbed tables
(Color difference on Low Base for illustration purposes only)

Specifications:

		OBT230 PatientMate®	OBT220 PatientMate® Jr.	OBT210 Standard	OBT100 C - Base
Overall Height	(hi/hi)	46" / 117 cm	46" / 117 cm	44" / 112 cm	42.5" / 108 cm
	(lo/lo)	30" / 76 cm	30" / 76 cm	29" / 74 cm	28.5" / 72 cm
Overall Width	(closed)	15" / 38 cm	15" / 38 cm	15" / 38 cm	18" / 46 cm
	(extended)	30" / 76 cm	30" / 76 cm	NA	NA
Food Tray Width		14.5" / 37 cm	14.5" / 37 cm	NA	NA
Overall Depth		44" / 112 cm	36" / 91 cm	36" / 91 cm	31" / 79 cm
Shipping Weight	(unpacked)	69 lbs / 31 kg	58 lbs / 26 kg	40 lbs / 18 kg	
	(packed)	79 lbs / 36 kg	67 lbs / 30 kg	48 lbs / 22 kg	
Standard Base		3" H x 19" W x 31" L			
Low Base			2" H x 19" W x 35" L		

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2007 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.

CTG164rb Rev. 2 3/09/07